

**NICKEL TITANIUM STEPPED FORCE ARCHWIRES**



Available in True Form - our most popular shape.

**It's Like Having Three Archwires in One!**

You treat your patients as individuals. Now you can apply the same philosophy to their teeth.

**Fewer wire changes required...**

The differentiated force values provide for levels of anterior/posterior correction previously only possible when using progressively heavier wires!

**Extend the intervals between appointments...**

Up to 10 weeks without loss of corrective potential of the archwires.

**ANTERIOR REGION – LIGHT FORCE**


ProFlex III™ delivers the lowest possible heat-activated force to the teeth with the smallest root structures. Wires are easily inserted, patients are more comfortable, and the risk of root resorption is greatly reduced.

**BICUSPID REGION – MEDIUM FORCE**

ProFlex III™ delivers progressively increased force to the bicuspid region for faster leveling and aligning versus conventional initial archwires.

**POSTERIOR REGION – HEAVY FORCE**

ProFlex III™ concentrates its greatest force in this region. Difficult molar rotations are corrected far more quickly, saving valuable chair time associated with multiple archwire changes.



10 Pk	Force Levels in Grams			Stepped-Force Archwire		
	Round	Anterior	Bicuspid	Posterior	Upper	Lower
.014"		30	45	60	W83U-014	W83L-014
.016"		50	60	90	W83U-016	W83L-016
.018"		60	75	100	W83U-018	W83L-018

Rectangular	Force Levels in Grams			Stepped-Force Archwire	
	Anterior	Bicuspid	Posterior	Upper	Lower
.016" x .016"	50	70	100	W83U-1616	W83L-1616
.016" x .022"	60	100	150	W83U-1622	W83L-1622
.017" x .025"	90	160	250	W83U-1725	W83L-1725
.018" x .018"	60	80	150	W83U-1818	W83L-1818
.018" x .025"	100	170	260	W83U-1825	W83L-1825
.019" x .025"	120	200	280	W83U-1925	W83L-1925
.020" x .020"	90	150	200	W83U-2020	W83L-2020

Wire Products